# SMTube®

## Strip Meniscometry Tube

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
Figure 1 : Summary of Subjects

Diagnosis	Subjects (eyes)	Gender *		Age Range (yrs)	Mean Age* (yrs)
		Men (eyes)	Women (eyes)		
DE	23 (45)	5 (10)	18 (35)	8-86	56.4 ± 17.0
Normal Controls	29 (51)	7 (11)	22 (40)	7-82	53.6 ± 18.4

Subjects with dry eye symptoms, presence of a vital staining score of >3points, and decreased tear break up time or tear quantity were diagnosed as having definite dry eye disease according to the 2006 Japanese Dry Eye Research Society Diagnostic Criteria<sup>4</sup>

\*No statistically significant differences were found in terms of male-female ratio ( $P = 1.0$  with Fisher's exact test,  $P = 0.90$  with the chi-square test with Yates' correction) and the age distribution ( $P = 0.61$  with Mann-Whitney's U test).

Figure 2 : Correlation among SMTube and other ocular surface examinations

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- PCT patented
- European Medical Device Directive 93/42/EEC, CE marked as Class I
- FDA device registered as Class I
- MHLW Japan registered as Class I medical device

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- 1. Screening of a number of subjects such as group health checkup. (Ishikawa 2016)**
- 2. Examining the amount of tears before and after a certain treatment and medication. (Ibrahim 2012, Kojima 2014, Okura 2015)**
- 3. Monitoring the time-series variation of the tear volume, which is repeated multiple times in relatively short time. (Kojima 2011, Hirayama 2013, Sano 2015)**
- 4. It is also suitable for tear volume evaluation in pediatrics, owing to its minimal invasiveness.**

# Références

## Original Papers of Strip Meniscometry

Dogru M, Ishida K, Matsumoto Y, et al. Strip meniscometry: a new and simple method of tear meniscus evaluation. *Invest Ophthalmol Vis Sci* 2006;47:1895-1901.

Ibrahim OM, Dogru M, Ward SK, et al. The efficacy, sensitivity, and specificity of strip meniscometry in conjunction with tear function tests in the assessment of tear meniscus. *Invest Ophthalmol Vis Sci* 2011;52:2194-2198.

## Publications involving SMTube product in their studies

Ishikawa S, Kato N, Okamoto Y, Takeuchi D. 健康診断におけるStrip Meniscometry を用いた涙液量検査の有用性: Utility of tear volume measurement using strip meniscometry at a mass medical checkup. *Annual Meeting of Japanese Society of Lacrimal Passage and Tear Dynamics*. Tokyo: Session Number RR-7-1; 2014. (in Japanese)

Shinzawa M, Dogru M, Miyasaka K, Shimazaki J, Sekiryu T. 前眼部OCT と strip meniscometry のドライアイ診断への応用: Application of Anterior-Segment OCT and strip meniscometry in the diagnosis of dry eye. *J Jpn Ophthalmol Soc* 2015;119(Suppl.):212. (in Japanese)

Sagara H. Quantitative Measurement of Tear Meniscus Volume during Visual Field Testing with Strip Meniscometry. *6th World Glaucoma Congress*. Hong Kong: Session Number P-T-057; 2015.

Sano K, Kawashima M, Ikeura K, Arita R, Tsubota K. Abdominal breathing increases tear secretion in healthy women. *Ocul Surf* 2015;13:82-87.

Okura M, Kawashima M, Katagiri M, Shirasawa T, Tsubota K. New Eye Cleansing Product Improves Makeup-Related Ocular Problems. *J Ophthalmol* 2015;2015:831628.

Ishikawa S, Kato N, Takeuchi M. The Efficacy of Strip Meniscometry and DEQS for Dry Eye Screening during Medical Health checkups. Poster Session in ARVO (Poster Board Number: A0061), 2016.

Shinzawa M, Dogru M, Miyasaka K, Shimazaki J, Sekiryu T. Application of CASIA SS-1000 Optical Coherence Tomography Tear Meniscus Imaging in Testing the Efficacy of New Strip Meniscometry in Dry Eye Diagnosis. *Eye Cont Lens* 2017.

Lee KW, Kim JY, Chin HS, Seo KY, Kim TI, Jung JW. Assessment of the Tear Meniscus by Strip Meniscometry and Keratograph in Patients With Dry Eye Disease According to the Presence of Meibomian Gland Dysfunction. *Cornea* 2017;36:189-195.

## Strip Meniscometry related papers

Dogru M, Ward SK, Wakamatsu T, et al. The effects of 2 week senofilcon-A silicone hydrogel contact lens daily wear on tear functions and ocular surface health status. *Cont Lens Anterior Eye* 2011;34:77-82.

Kojima T, Matsumoto Y, Ibrahim OM, et al. Effect of controlled adverse chamber environment exposure on tear functions in silicon hydrogel and hydrogel soft contact lens wearers. *Invest Ophthalmol Vis Sci* 2011;52:8811-8817.

Ibrahim OM, Dogru M, Kojima T, et al. OCT assessment of tear meniscus after punctal occlusion in dry eye disease. *Optom Vis Sci* 2012;89:E770-E776.

Hirayama M, Murat D, Liu Y, Kojima T, Kawakita T, Tsubota K. Efficacy of a novel moist cool air device in office workers with dry eye disease. *Acta Ophthalmol* 2013;91:756-762.

Ibrahim OM, Dogru M, Kawashima S, Wakamatsu TH, Tsubota K, Fujishima H. Visante optical coherence tomography and tear function test evaluation of cholinergic treatment response in patients with sjogren syndrome. *Cornea* 2013;32:653-657.

Kojima T, Matsumoto Y, Ibrahim OM, Wakamatsu TH, Dogru M, Tsubota K. Evaluation of a thermosensitive atelocollagen punctal plug treatment for dry eye disease. *Am J Ophthalmol* 2014;157:311-317.

## PRODUCT INQUIRIES : TOKYO OFFICE

6-24-14 Umeda Adachi-ku Tokyo 123-0851 Japan  
PHONE : +81 3 3880 0455 FAX : +81 3 3848 1825  
URL : <http://www.echo-mf.jp/en.html>  
EMAIL : [contact-english@echo-mf.jp](mailto:contact-english@echo-mf.jp)

## MANUFACTURER : SHIRAKAWA FACTORY

1-2 Shinnatsunashi Shirakawa-shi Fukushima-ken 961-0024 Japan

- ISO 13485:2003, ENISO13485:2012 certified
- Medical devices manufacturer licensed by MHLW Japan
- MHLW third-class medical devices marketing authorization holder
- FDA registration