

## Coaxial Disposable Irrigation Aspiration System 20G, 21G, 22G




### Product Features:

- Coaxial System Disposable 20G & 22G, 21G optional
- Lightweight irrigation / aspiration Handpiece
- Prolonged sleeve (plus 1.5mm)
- all G: 45° tip smooth abrasive (rough) or polished
- 20G, 21G: 65° tip smooth abrasive (rough) or polished
- Compatible with your existing Phaco or I/A System
- Swiss Design / Swiss Made


### Details:




prolonged tip



45° tip



65° tip



./.

**CX2045-P**

Coaxial System Disposable 20G  
Asp. port 0.35mm, angled 45°, POLISHED

**CX2065-P**

Coaxial System Disposable 20G  
Asp. port 0.35mm, angled 65°, POLISHED

**CX2145-P**

Coaxial System Disposable 21G  
Asp. port 0.35mm, angled 45°, POLISHED

**CX2165-P**

Coaxial System Disposable 21G  
Asp. port 0.35mm, angled 65°, POLISHED

**CX2245-P**

Coaxial Sytem Disposable 22G  
Asp. port 0.35mm, angled 45°, POLISHED

**CX2045-R**

Coaxial System Disposable 20G  
Asp. Port 0.35mm, angled 45°, SMOOTH  
ABRASIVE (ROUGH)

**CX2065-R**

Coaxial System Disposable 20G  
Asp. port 0.35mm, angled 65°, SMOOTH  
ABRASIVE (ROUGH)

**CX2145-R**

Coaxial System Disposable 21G  
Asp. port 0.35mm, angled 45°, SMOOTH  
ABRASIVE (ROUGH)

**CX2165-R**

Coaxial System Disposable 21G  
Asp. port 0.35mm, angled 65°, SMOOTH  
ABRASIVE (ROUGH)

**CX2245-R**

Coaxial Sytem Disposable 22G  
Asp. port 0.35mm, angled 45°, SMOOTH  
ABRASIVE (ROUGH)